



> Product Group Laparoscopy

- A.M.I.-dock Instruments
- A.M.I. Disposable Instruments
- A.M.I. Knotpusher Premium
- A.M.I. EasyInstruments
- A.M.I. OptiClean Device
- A.M.I. SituClean Device
- A.M.I. EasyTac System
- A.M.I. Hernia Mesh

> We are at your service:

A.M.I. Ltd.

Im Letten 1 . 6800 Feldkirch . Austria




t +43 (0) 5522 90505-0

f +43 (0) 5522 90505-4006

e customerservice@ami.at

i www.ami.at

A.M.I.[®]-dock Instruments

Order Code	
A05 491	A.M.I.-dock Applicator Standard 5 mm
A10 111	A.M.I.-dock Applicator Standard 10 mm
	<p>Working length: 510 mm The universal instrument in minimal invasive surgery Reusable stainless steel instrument Non sterile, steam autoclavable. 1 unit</p>
A05 492	A.M.I.-dock Applicator Short 5 mm
A10 112	A.M.I.-dock Applicator Short 10 mm
	<p>Working length: approx. 390 mm Reusable stainless steel instrument for paediatric surgery. Non sterile, steam autoclavable. 1 unit</p>
A05 493	A.M.I.-dock Applicator Curved 5 mm
	<p>Working length: 510 mm Reusable stainless steel instrument for thoracic surgery. Non sterile, steam autoclavable. 1 unit</p>
A05 891	Rubber ring seal 5 mm
A10 892	Rubber ring seal 10 mm
	<p>To be used in combination with the A.M.I.-dock Applicator Standard or Short. Non sterile, steam autoclavable. 2 units / box</p>

A.M.I.[®]-dock Disposable Loading Units

Order Code	Product
A10 211	A.M.I.-dock TissueBag DLU For safe removal of specimen (max. 65 mm in diameter) during laparoscopic procedures. Disposable Loading Unit to be used with the A.M.I.-dock Applicator 10 mm, delivered sterile. 10 units / box
A10 212	A.M.I.-dock TissueBag Large DLU For safe removal of specimen (max. 100 mm in diameter) during laparoscopic procedures. Disposable Loading Unit to be used with the A.M.I.-dock Applicator 10 mm, delivered sterile. 10 units / box
A05 251	A.M.I.-dock BluntDissect DLU 5 mm
A10131	A.M.I.-dock BluntDissect DLU 10 mm For blunt dissection during minimal invasive procedures. Disposable Loading Unit to be used with the A.M.I.-dock Applicator, delivered sterile. 20 units / box
A05 451	A.M.I.-dock Easy Loop DLU 5 mm Tissue bundles and vessels can be ligated under enhanced control due to loop fixation. Disposable Loading Unit to be used with the A.M.I.-dock Applicator 5 mm, delivered sterile. 20 units / box
A05 321	A.M.I.-dock Pleuroscratch DLU 5 mm Disposable Loading Unit to be used with the A.M.I.-dock Applicator 5 mm, delivered sterile. Easy and controlled traumatisation of parietal pleura. 20 units / box

A.M.I.® Disposable Instruments

Order Code	Product
TBI 100	A.M.I. TissueBag Disposable Instrument Fully disposable Instrument for the safe removal of specimen (max 65 mm in diameter) The instrument is considered for a 10 mm trocar. Delivered sterile. 10 units / box
AFC 100	A.M.I. FasciaCutter Disposable Instrument Instrument with a short handle for safe enlargement of trocar incisions. Delivered sterile 20 untis / box

A.M.I.® Laparoscopic Knotpusher

Order Code	Product
APK 001	A.M.I. Knotpusher Premium Reusable Stainless Steel Instrument for laparoscopic surgery. Length of shaft: 400 mm Length of instrument: 500 mm Non sterile, steam autoclavable 1 unit

A.M.I.® EasyShears Instrument






Order Code	Product
ESI1031	A.M.I. Disposable Tips  To be attached to the A.M.I. EasyShears Instrument Total length of the curved blade: 16 mm Delivered sterile. 10 units / box
ESI1051	A.M.I. Mini Disposable Tips  To be attached to the A.M.I. EasyShears Instrument Total length of the curved blade: 10 mm Delivered sterile. 10 units / box
ESI1020	A.M.I. Reusable Tips  To be attached to the A.M.I. EasyShears Instrument Total length of the curved blade: 16 mm Delivered non sterile. Steam autoclavable. 5 units / box
ESI1040	A.M.I. Mini Reusable Tips  To be attached to the A.M.I. EasyShears Instrument Total length of the curved blade: 10 mm Delivered non sterile. Steam autoclavable. 5 units / box

Please take note:

These tips are considered to be attached to the A.M.I. EasyShears ESI1010.

They do not match to the A.M.I. Easy Instrument Premium ESI1110.

A.M.I.® Easy Instrument Premium

Order Code	Product
<p data-bbox="185 327 280 349">ESI1110</p> 	<p data-bbox="451 327 839 349">A.M.I. Easy Instrument Premium</p> <p data-bbox="451 387 1422 533">Reusable Instrument 5 mm, shaft 330 mm long, fully detachable, 360° shaft rotation, exchangeable instrument tips. Delivered non sterile, steam autoclavable. Can be combined with all kind of _ A.M.I. EasyShear Tips Premium, as well as with the _ A.M.I. EasyDissect Tips Premium</p> <p data-bbox="451 562 520 584">1 unit</p>
<p data-bbox="185 680 280 703">ESI1131</p> 	<p data-bbox="451 680 967 703">A.M.I. EasyShear Disposable Tips Premium</p> <p data-bbox="451 741 1046 819">Attachable to the A.M.I. Easy Instrument Premium Total length of the curved blade: 16 mm Delivered sterile.</p> <p data-bbox="451 853 616 875">10 units / box</p>
<p data-bbox="185 945 280 967">ESI1151</p> 	<p data-bbox="451 945 1026 967">A.M.I. EasyShear Mini Disposable Tips Premium</p> <p data-bbox="451 1005 1046 1084">Attachable to the A.M.I. Easy Instrument Premium Total length of the curved blade: 10 mm Delivered sterile.</p> <p data-bbox="451 1117 616 1140">10 units / box</p>
<p data-bbox="185 1209 280 1232">ESI1120</p> 	<p data-bbox="451 1209 943 1232">A.M.I. EasyShear Reusable Tips Premium</p> <p data-bbox="451 1270 1046 1382">Attachable to the A.M.I. Easy Instrument Premium Total length of the curved blade: 16 mm Delivered non sterile. Steam autoclavable.</p> <p data-bbox="451 1415 600 1438">5 units / box</p>
<p data-bbox="185 1505 280 1527">ESI1140</p> 	<p data-bbox="451 1505 1002 1527">A.M.I. EasyShear Mini Reusable Tips Premium</p> <p data-bbox="451 1565 1046 1677">Attachable to the A.M.I. Easy Instrument Premium Total length of the curved blade: 10 mm Delivered non sterile. Steam autoclavable.</p> <p data-bbox="451 1711 600 1733">5 units / box</p>

A.M.I.® Easy Instrument Premium

Order Code	Product
ESI1110	<p>A.M.I. Easy Instrument Premium</p> <p>Reusable Instrument 5 mm, shaft 330 mm long, fully detachable, 360° shaft rotation, exchangeable instrument tips. Delivered non sterile, steam autoclavable. Can be combined with all kind of _ A.M.I. EasyShear Tips Premium, as well as with the _ A.M.I. EasyDissect Tips Premium 1 unit</p>
EDI1031	<p>A.M.I. EasyDissect Disposable Tips Premium</p> <p>Single use loading unit for dissection to be attached to the A.M.I. Easy Instrument Premium. Delivered sterile. 10 units / box</p>
EDI1020	<p>A.M.I. EasyDissect Reusable Tips Premium</p> <p>Multiple use loading unit for dissection to be attached to the A.M.I. Easy Instrument Premium. Delivered non sterile. Steam autoclavable. 5 units / box</p>

Order Code	Product
ECI1010	<p>A.M.I. EasyClinch Instrument Premium</p> <p>Reusable Instrument 5 mm, shaft 330 mm long, fully detachable, 360° shaft rotation, exchangeable instrument tips. Can be combined with A.M.I. EasyClinch Tips Premium. Delivered non sterile, steam autoclavable. 1 unit</p>
ECI1020	<p>A.M.I. EasyClinch Reusable Tips Premium</p> <p>Atraumatic grasper. Multiple use loading unit to be attached to the A.M.I. EasyClinch Instrument Premium. Delivered non sterile. Steam autoclavable. 5 units / box</p>

A.M.I.® OptiClean System

Order Code	Product
AOC 400	<p>A.M.I. OptiClean Device 305</p> <p>Disposable device to keep the lens of any scope clean during endoscopic procedures, without the need to remove the scope from the surgical site</p> <p>System consists of:</p> <ul style="list-style-type: none"> - A.M.I. OptiClean Snap on Profile, Length 305 mm - Tubing - Irrigation Syringe - Suction tube puncture needle <p>6 units / box</p>
AOC 402	<p>A.M.I. OptiClean Device 285</p> <p>Disposable device to keep the lens of any scope clean during endoscopic procedures, without the need to remove the scope from the surgical site</p> <p>System consists of:</p> <ul style="list-style-type: none"> - A.M.I. OptiClean Snap on Profile, Length 285 mm - Tubing - Irrigation Syringe - Suction tube puncture needle <p>6 units / box</p>

A.M.I.® SituClean Device

Order Code	Product
ASC 420	<p>A.M.I. SituClean Device</p> <p>Suction / Irrigation Device</p> <p>For rinsing the surgical site and removal of small volumes of blood and body liquids during laparoscopic or arthroscopic procedures</p> <p>Sterile packed</p> <p>3 units / box</p>

A.M.I.® EasyTac System

Order Code	Product
ETI 100	<p>A.M.I. EasyTac Instrument</p> <p>Length of shaft: 297 mm Reusable stainless steel instrument 5 mm, fully detachable. Use in connection with the A.M.I. EasyTac DLU for mesh fixation in hernia repair operations. Non sterile, steam autoclaveable. 1 unit</p>
ETI 110	<p>A.M.I. EasyTac Instrument Short</p> <p>Length of shaft: 97 mm Reusable stainless steel instrument 5 mm, fully detachable. Use in connection with the A.M.I. EasyTac DLU for mesh fixation in hernia repair operations. Non sterile, steam autoclaveable. 1 unit</p>
ETD 121	<p>A.M.I. EasyTac DLU</p> <p>Disposable Loading Units, 5 mm Including 5 fixation elements Use in connection with the A.M.I. EasyTac Instrument and the A.M.I. EasyTac Instrument Short Sterile packed. 5 units / box</p>
ETI 121	<p>A.M.I. EasyTac Disposable Instrument Premium</p> <p>Length of shaft: 383 mm Fully disposable instrument 5 mm, including 12 fixation elements Use for mesh fixation in hernia repair operations Sterile packed. 1 unit</p>
ETI 131	<p>A.M.I. EasyTac Disposable Instrument Short Premium</p> <p>Length of shaft: 183 mm Fully disposable instrument 5 mm, including 12 fixation elements Use for mesh fixation in hernia repair operations Sterile packed. 1 unit</p>

A.M.I.® EasyTac System

Order Code	Product
ETI 000	A.M.I. EasyTac DEMO Demo System consists: - A.M.I. EasyTac Disposable Instrument Short Length of shaft: 183 mm, including 5 fixation elements, 1 unit - Mesh Substitute, 30 x 30 mm - Foam Block, 50 x 50 mm NON STERILE - FOR DEMONSTRATION ONLY

A.M.I.® Hernia Mesh

Order Code	Product
AHM 1011	A.M.I. Hernia Mesh Size: 150 x 150 mm Monofilament, ultralight Polypropylene, 19 g/m ² 5 units / box
AHM 1012	A.M.I. Hernia Mesh Size: 300 x 300 mm Monofilament, ultralight Polypropylene, 19 g/m ² 5 units / box